

# POLISH MANUFACTURER OF SPECIAL NEEDS STROLLERS

www.mywam.eu



NOTE: The information in this leaflet is provided for convenience and education purposes only and therefore: a) such information does not constitute a commercial proposal (quotation) within the meaning of commercial law, b) such information cannot be the basis for taking decisions about taking up treatment, withdrawal from treatment or other medical procedures, c) the information should be interpreted by a medical professional. d) the products may be modified by the manufacturer without notice.

This is a medical device.

Use in accordance with the instruction manual or label.





### POLISH MANUFACTURER OF SPECIAL NEEDS STROLLERS

C € MD MEDICAL DEVICE - CLASS I

The MOUSE special needs stroller is designed for the smallest children who are unable to walk or sit in a correct position as a result of injuries, diseases or motor dysfunctions. The stroller is suitable for children at a height of 80-110 cm and a weight of 30 kg.

Thanks to the backrest tilt and recline feature as well as adjustable torso and head pads you may easily adjust MOUSE to the required position and ensure the user with optimal protection and stability.

The shock-absorbing suspension on all four, large, gel-filled wheels together with the tilt-foot facilitate travel on uneven terrain. Modern and stylish design is an additional advantage of this stroller.

#### **INDICATIONS**

paralysis and paresis of lower limbs of various origins, paralysis and paresis of the trunk, weakening of postural muscles responsible for maintaining the correct body position, cerebral palsy, myelomening ocele, muscular dystrophies.

#### CONTRAINDICATIONS:

In case of any contraindications to sitting position &/or passive verticalization stroller adjustments should be consulted with physiotherapist or physician-in-charge.

# colour options.



FRAME
spray painted / glossy











# features. safety and comfort

**TILT IN SPACE SEAT** 

simple mechanism



8°-23°

4 positions

**RECLINING BACKREST** wide range of adjustments



16°- 88°

4 positions

### ADJUSTABLE AND ANGLE TILT FOOTREST

simple and quick adjustment



84°- 173°

6 positions

#### **SEAT MINIMIZING KIT**

reducing insert



The SEAT MINIMIZING KIT in MOUSE is a comfortable, one-peace insert that allows you to reduce both width and depth of the seat.

#### STABILIZATION SYSTEM

adjustable pads



The head and torso support consists of two pairs of independent pads that are height, width and angle adjustable. Its wide range of adjustments and simplicity gives you a fast and effective adaptation to the user's individual needs.

#### LYING POSITION

wide range of backrest recline



With a fully adjustable seat and footrest you may quickly alter the stroller into a cozy bed perfect for a short rest or a nap.

#### **REVERSIBLE SEAT**

arrange your travel easily



The design of the stroller allows you to easily rotate the seat opposite to the driving direction. Thanks to this feature you can keep an eye on your child all the time

#### **FOLDING LIGHTWEIGHT FRAME**

ease of storage, transportation and dismantling



The construction of main elements allows you to dismantle Mouse into 3 separate units (frame, seat and wheels).

#### **DRIVING SUPPORT**

comfort for the bumpiest journeys



The stroller is equipped with a shock-absorbing suspension on all four wheels and an adjustable rear suspension in the rear wheels, swivel front castor wheels with option to lock, central brake and atilt-foot.

# functional accessories. in practice



## accessories. wide range of equipment



### parameters.

#### TECHNICAL DATA OF A STROLLER READY TO DRIVE



seat depth 22\*/26 cm

backrest height 52 cm

footrest adjustable within 14,5 - 23 cm

**6** total weight 17,2 - 23,2 kg

total lenght 95 - 124 cm

total height 122 - 134 cm

total width 60 cm

dimension after folding: frame L 84 x W 58 x H 25 cm / seat L 72 x W 42 x H 39 cm

#### **USER PARAMETERS**

🙀 growth 80 - 110 cm

maximum weight 30 kg